



MULTI-RESIDUE PESTICIDE ANALYSIS TECHNIQUES, AN OVERVIEW

PMN Rajesh

Instruments and workflow



SAMPLE EXTRACTION WORKFLOW



Waters Product line segmentation Food Testing

Sample
preparation



GC-MS,APGC
Systems

Application
Managers
& Informatics



TargetLynx™



UNIFI

Progenesis® Q1
Nonlinear
dynamics



MS
Detectors
Tandem
Quad



HRMS - Q TOF



Why mass spectrometry required

Spectrometry - LC/UV/FL, HPTLC

- Insensitive, non- specific
- Rigorous sample prep and derivatization often required
- Difficulty to analyze multiple analytes in a single run


Immunoassay - RIA, EIA (screening only)

- Not specific, issues with cross contamination, false positive
- Expensive and only available for a few compound classes
- No way to analyze multiple analytes in a single run

LC/MS/MS

- Highly selective and accurate
- Reduced sample preparation
- No issue with false positives
- Provides both qualitative and quantitative analysis in a single run
- Meet all requirements for modern residue analysis
- Technique recommended by all food safety guidelines

What are the challenges?



What impact does this have on the laboratory?

- Demand for Flexible methods
 - Many analytes
 - Many matrices
 - Simple set-up / modification
- Method performance
 - High sensitivity
 - Quantitative accuracy
- Data management
 - Automated processing
 - Reporting
 - Archive and back-up



Harmonisation of analytical quality control and validation

- **SANTE/12682/2019** EU Guidelines on Analytical Quality Control and Validation Procedures for Pesticide Residues Analysis in Food and Feed
- Provides detailed advice on **best laboratory practice** from sample receipt though to reporting of the result(s)
- Includes **identification** criteria (NOT identification points)
 - Allows the use of accurate mass
- Intended as **guidance** for official control purposes
- **Regularly reviewed and updated**

Selection of representative commodities for validation

- High water content
 - Apples, carrots, lettuce, peaches, tomatoes
- High acidity (and high water content)
 - Berries, citrus, grapes, pineapples
- High sugar (and low water content)
 - Dried fruits, jams
- High starch and/or protein content (low water and fat content)
 - Cereals, dried beans pulses
- High oil (and very low water content)
 - Tree nuts, oil seeds
- High oil (and intermediate water content)
 - Avocados, olives
- Unique/difficult
 - Cocoa, coffee, hops, spices, teas
- Meat and seafood
 - Beef, chicken muscle, liver, kidney, cod, salmon, shrimp
- Milk and milk products
 - Milk, cheese
- Eggs
- Fat from food of animal origin
 - Lard, butter, cod liver oil

XEVO Tandem Quadrupole

XEVO is

- ✓ Robust, reliable and accessible high performance
- ✓ Innovative, information rich analysis
- ✓ Versatility with engineered simplicity

● **Xevo TQD**



**Accessible,
reliable and
proven**

● **Xevo
TQ-S micro**



**Sensitive, reliable
and compact**

● **Xevo TQ-S**



**High sensitivity and
reliability**

● **Xevo TQ-XS**




**Ultimate sensitivity
and robustness**


- **Increasing Instrument Uptime**
- **Reducing Sample Turnaround Time**
- **Easy Method Transfer**
- **Understanding Sample Complexity**

ESCi Multi-mode Ionisation

- ESI and APCI scans in the same acquisition
- Allows wider range of compounds to be ionised




ESCi mode




Xevo TQ-S micro Innovative Technologies


■ **STEP WAVE™** ion guide



■ **SPEED** with high quality data



■ **RADAR™** and PICs



■ Quantitative workflow updates
– Standard Addition
– Increased no. compounds



Informatics solutions

INTELLIⁱSTART™

QUANPEDIA™

QCMonitor

TargetLynx™ XS

Trendplot

- Input the compound name
- Calculate potential precursor masses
- Select Multiply Charged Parents
- Input the multiply charged precursor mass(es) from MW calculator
- Select advanced mode
- Limit low mass fragments
- Click Start

IntelliStart Setup Parameters

Sample Tune and Develop Method

Compound Details

	Compound Name	Parent Mass	Ion Mode
<input checked="" type="checkbox"/>	Insulin	1434	ES+ ▼
<input checked="" type="checkbox"/>	Insulin 5+	1148	ES+ ▼
<input checked="" type="checkbox"/>	Insulin 6+	956	ES+ ▼
<input type="checkbox"/>			ES+ ▼

☒ Multiply Charged Parents

Method Details

☒ Create New Sample Tune ☐ Load Existing Sample Tune

Sample Tune Name: Insulin.IPR

Develop SIR method: SIR_Insulin

Develop MRM method: MRM_Insulin

Save Report As: Report_Insulin

☐ Invoke Manual Optimisation

☒ Append to existing methods

☐ Export To LC/MS System Check

☒ Print Report

Optimization Ranges

Cone Voltage: Default (2 - 100) V

Collision Energy: Peptide (2 - 50) eV

Daughter ion settings

Number of MRM transitions per compound: 5


Lowest Fragment Ion Mass: 200.0 Da

Fluidics


Flow Path: Combined

Reservoir: A


Sample Solution



TargetLynx






Rapid MS to MRM Switching




Knowing more about your samples RADAR

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Information rich data


Precursor Ion Confirmation Scan



“Dilute and shoot” 10 ppb pesticides



ROTENONE





Sensitivity and Matrix effect



Methamidophos @ 10 ppb




10x dilution



Quechers
extract
10 %
CH₃CN



TargetLynx XS Standard Addition



SAMPLE EXTRACTION

Multi Residue Analysis for Food Safety



Typical challenges include


- ? Sample throughput: >1000 pesticides commercially available foods and feeds
- ? Regulatory requirements: Approx 740 compounds regulated concentrations
- ? Accessibility of analysis: LC and GC amenable pesticides
- ? Nontargeted analysis: Metabolites or unknown contaminants

PASS ☒
FAIL ☐
SANTE 11813/2017

DisQuE Products




Basic QuEChERS Procedure



Extraction/Partitioning

CEN 15662 Method



- Transfer 10 g of sample to a 50 mL extraction tube
- Add 10 mL acetonitrile. Shake for 1 minute
- Add contents of DisQuE pouch
- Shake vigorously for 1 minute
- Centrifuge, take aliquot of top layer for cleanup by dSPE or other methods

QuEChERS Extraction: 50 mL Tube

■ Water is IMPORTANT!!


- Produces liquid/liquid partition cleanup after DisQuE salts are added
 - Removes unwanted polar matrix components (sugars, salts)
 - Makes acetonitrile extraction more effective
- If water not present in matrix, then add it

Acetonitrile Layer (**ANALYTES**)

Remaining sample solids

Aqueous Layer (**Saturated buffer salts
+ ionic/polar analytes**)

Undissolved buffer salts



Quechers

QuEChERS extraction



QuEChERS
extract



QuEChERS extract with
dSPE cleanup/10X



QuEChERS
extract 20X



New DisQuE Cleanup Tubes

- Available Now

Tube Size	Contents	Description
2 mL	150 mg MgSO ₄ , 25 mg PSA, 25 mg C18, 7 mg GCB	General purpose cleanup
15 mL	1200 mg MgSO ₄ , 400 mg PSA	AOAC Method
15 mL	1200 mg MgSO ₄ , 400 mg PSA, 400 mg C18	AOAC Method, fatty matrix
15 mL	1200 mg MgSO ₄ , 400 mg PSA, 400 mg C18, 400 mg GCB	AOAC Method, fatty & pigment matrix
2 mL	150 mg MgSO ₄ , 50 mg C18	Acidic pesticides
2 mL	150 mg MgSO ₄ , 25 mg PSA, 2.5 mg GCB	CEN Method, pigment matrix
15 mL	900 mg MgSO ₄ , 300 mg PSA	High fatty acids samples
15 mL	900 mg MgSO ₄ , 300 mg PSA, 300 mg C18	Fatty produce/cereal
15 mL	900 mg MgSO ₄ , 450 mg PSA, 300 mg C18, 50 mg GCB	Fatty and pigment matrix, Tea
15 mL	750 mg MgSO ₄ , 250 mg PSA, 150 mg C18, 150 mg Al-N	Multi-residue mycotoxins
2 mL	150 mg MgSO ₄ , 50 mg PSA, 30 mg C18, 30 mg Al-N	Multi-residue mycotoxins

✓ Consider Oasis PRiME HLB pass-thru cleanup

Other Cleanup Options

- by Cartridge SPE

- Sep-Pak Vac 6 cc (500mg) PSA/Carbon cartridges
 - Used for APGC-MS cleanup of spices
- Oasis PRiME HLB for **pass-thru** cleanup
 - As effective as dSPE for removal of
 - Fats
 - Phospholipids
 - Chlorophyll
- Acidic pesticides – PSA cleanup is not suitable
 - Diluted with mobile phase and analyzed directly
 - If cleanup is required, e.g. tea
 - Use Oasis MAX 3 cc 60 mg for acidic pesticides cleanup in tea



SPE Procedure

Step 1

- Condition HLB and MCX the cartridges with 1.5mL of methanol and 1.5mL of water

Step 2

- 10mL of centrifuged wastewater is passed through the cartridges

Step 3

- Wash HLB with 3mL of water
- Wash MCX with 3mL of 2% Formic acid in water

Step 4

- Elute with HLB with 2mL methanol
- Elute with MCX with 2mL methanol with 2% NH₄OH

Step 4

- Pool all elutions and evaporate to dryness and reconstitute with 1mL of Mobile phase

Experimental Overview



Sample Preparation

- CEN QuEChERS, tube 1

Liquid Chromatography Conditions

- ACQUITY UPLC H-Class
- Column: ACQUITY UPLC BEH C₁₈ 100 mm x 2.1 mm, 1.7 μ m
- Column temperature: 45 ° C
- MP A: 10 mM ammonium acetate (pH 5) in water
- MP B: 10 mM ammonium acetate (pH 5) in acetonitrile
- Flow rate: 450 μ l
- Injection volume: 10 μ l

MS Conditions

- Ionisation: ES \pm
- Capillary: 1 kV (+) and 0.5 kV (-)
- Source Temperature : 150 °C
- Desolvation Temperature : 1000 °C
- Cone Gas flow: 120 l.hr⁻¹
- Desolvation Gas flow: 1000 l.hr⁻¹

Chromatogram of 400 pesticides



Waters
THE SCIENCE OF WHAT'S POSSIBLE.™






Figure 1. Chromatogram showing all 402 pesticide residues in one 10 minute run in injection solvent.



200 Pesticides with two transitions in 12 min on Xevo TQ-S micro




Chromatogram of Chilli extracted sample for five pesticides in dual mode




Chromatogram for Aflatoxin B1 (0.015ng/mL)



Chromatograms of Tetracycline's, Quinolone antibiotics and Sulfonamides in Single run



Nitro furan parent compounds



POLAR PESTICIDES

Sample Preparation for Polar pesticides

Approaching multi residue extraction with QuPPe method

Weigh homogenized sample (10g) into centrifuge tube
(adjust for water content)


← spike commodity with labelled IS or reference material

↓
Add methanol (10 ml) containing 1 % formic acid

↓
Vortex thoroughly for 1 minutes

↓
Centrifuge at 5000 rpm for 5 minutes

↓
Filter supernatant (0.45 μ m, PVDF, filter) into a plastic vial



**Quick Method for the Analysis of Numerous
Highly Polar Pesticides in Foods of Plant Origin via LC-MS/MS
Involving Simultaneous Extraction with Methanol (QuPPe-Method)
I. Food of Plant Origin (QuPPe-PO-Method)**

Version 10 (09.01.2019, Document History, see page 73)


Authors: M. Anastasiades; D. I. Kolberg; E. Eichhorn; A. Benkenstein; A.-K. Wachtler; S. Zechmann;
D. Mack; C. Wildgrube; A. Barth; I. Sigalov; S. Görlich; D. Dörk; G. Cerchia

EU Reference Laboratory for pesticides requiring Single Residue Methods (EURL-SRM)
Address: CVUA Stuttgart, Schaflandstr. 3/2, DE-70736 Fellbach, Germany
Web: www.eurl-pesticides.eu,
E-Mail: EURL@cvuas.bwl.de

Note: Changes from V9.3 to V10 are highlighted in yellow




Waters' polar pesticide journey simplified




Torus Technology

- HILIC type: DEA
- ✓ Retention
- ✓ Excellent separation of critical pairs
- ✓ Retention time stable
- ✓ Perchlorate compatible

Separation




Retention time stability Glyphosate



- HILIC type (polyvinyl alcohol)
- ✓ Retention
- ✓ Acceptable separation of critical pairs
- ✗ Retention time unstable
- ✗ Perchlorate incompatible (5 min wide peak)


- Mixed mode
- ✗ Retention
- ✗ Separation of critical pairs

Retention



RSDs > 20% in different commodities

Chromatogram for Anionic polar pesticides analysis




Atmospheric Pressure Gas Chromatography




- APGC is an inlet option of the Xevo universal source, readily interchangeable with UPLC, *etc.*
- Complementary technique to GC-EI
- Operates
 - At atmospheric pressure, thus allowing higher gas flows to be applied
 - By APCI like ionisation using a corona discharge pin, thus providing 2 types ionisation

APGC – How It Works





Charge Transfer




“Dry” source conditions

Favoured by relatively **non-polar compounds**




Increased precursor ion providing confidence



Softer ionisation, improving selectivity: GC-EI-MS vs GC-APCI-MS/MS


Fipronil

Fipronil	EI	APGC
Quan	367 > 255	436 > 368
Qual	367 > 213	369 > 215



Fipronil sulfone

Sulfone	EI	APGC
Quan	383 > 335	451 > 415
Qual	383 > 255	451 > 282



Increased selectivity improves sensitivity
10 times dilution of 0.005 mg/kg

APGC provides softer ionisation and thus more confidence with molecular ion [M+H]⁺ information

Experimental:


Xevo G2 XS QToF Parameters

GC-MS Ionisation

Atmospheric pressure (API)

Corona current (μA)	3.0
Sampling Cone	20.0
Source Temperature	120° C
Source Offset	80
Cone Gas Flow	175 L/Hr
Auxiliary Gas Flow	50 L/Hr
Acquisition range	50-1200 m/z
Scan time	0.25 sec
Lockmass	Siloxane Bleed

Protonation



LC-MS Ionisation

Electrospray (ESI)

Capillary (kV)	0.8
Sampling Cone	20.0
Source Temperature	120° C
Source Offset	80
Desolvation Temperature	550° C
Cone Gas Flow	50 L/Hr
Desolvation Gas Flow	1000 L/Hr
Acquisition range	50-1200 m/z
Scan time	0.25 sec
Lockmass	LeuEnk (556.2771m/z)





Targeted Screening



Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Component Summary												
Component name	Identification status	Formula	Neutral mass (Da)	Observed neutral mass (Da)	Observed m/z	Mass error (mDa)	Mass error (ppm)	Observed RT (min)	Response	Adducts	Isotope Match	Mz RMS PPM
148 Hexazinone	Identified	C12H20N4O2	252.15863	252.1575	253.1647	-1.2	-4.6	10.72	304933	+H		4.54
149 Iodofenfos	Identified	C8H8Cl2IO3PS	411.83535	411.8352	412.8424	-0.2	-0.5	10.12	352525	+H		3.93
150 Iprodione	Identified	C13H13Cl2N3O3	329.03340	329.0341	330.0413	0.7	2.0	10.10	2360	+H		5.93
151 Isazophos	Identified	C9H17ClN3O3PS	313.04168	313.0401	313.0396	-1.5	-4.9	10.91	9993	-e		



Accurate Mass HRMS: Suspect to unknown screening

Are these compounds in my sample?

Screening

How much is in my sample?

What else is in my sample?

What is the difference between my sample and another one?

Comparison



Elucidation Toolset

Discovery

Halogen
Match

Common
Fragment
Match

Neutral
Loss

Mass Defect
Filter

UNIFI™

Identification

Library Search
(UNIFI
Scientific
Library)


Chemspider

Isotope
Match

Isotope
Model

Elemental
Composition

Qualitative Screen of Spinach: Binary Compare



Unknown Screening: Halogen Match Spinach

The screenshot shows the Waters software interface. The top navigation bar includes 'Review', 'Investigate', and 'Report' tabs. Below this is a 'Review Results' section with a search bar and various tool icons like 'Limit faults', 'Process', 'Edit', 'Tools', and 'File'. The main window displays a 'Component Summary' table with 8 rows of candidate masses. The left sidebar shows a 'Workflow' tree with 'Summary', 'Binary Compare 1_8', and 'Halogen Match Results' (highlighted). The 'Halogen Match Results' section shows a list of compounds: Ametoctradin, Bixafen, Penflufen, Pyriofenone, and Valifenalate.

Component name	Formula	Identification status	Expected RT (min)	Retention Time Error (min)	Mass error (ppm)	Response	Adducts
1 Candidate Mass 366.1108		None				951856	
2 Candidate Mass 399.1683		None				252365	
3 Candidate Mass 371.0622		None				136234	
4 Candidate Mass 414.0385		None				128171	
5 Candidate Mass 388.0922		None				109660	
6 Candidate Mass 421.1502		None				102809	
7 Candidate Mass 371.1006		None				81565	
8 Candidate Mass 485.1131		None				57547	

Spinach
(1.0 g/mL)
Ametoctradin
Bixafen
Penflufen
Pyriofenone
Valifenalate

Incurred/
comment
unk
unk
unk
unk
unk

Formula
 $C_{15}H_{25}N_5$
 $C_{18}H_{12}Cl_2F_3N_3O$
 $C_{18}H_{24}FN_3O$
 $C_{18}H_{20}ClNO_5$
 $C_{19}H_{27}ClN_2O_5$

Binary Compare and Elucidation
Binary Compare and Elucidation
Binary Compare and Elucidation
Binary Compare and Elucidation
Binary Compare and Elucidation

thank
you!

