

# Performance of ELISA-based Analytical Methods Supporting Compliance with Gluten-Free Regulations



# Gluten analysis and gluten-free products


# Gluten Analysis in Codex Alimentarius



## STANDARD FOR FOODS FOR SPECIAL DIETARY USE FOR PERSONS INTOLERANT TO GLUTEN

CODEX STAN 118-1979

Adopted in 1979. Amendment: 1983 and 2015. Revision: 2008.



## RECOMMENDED METHODS OF ANALYSIS AND SAMPLING

CXS 234-1999<sup>1</sup>

Adopted in 1999.

<sup>1</sup> The most updated version of the method should be used, in application of ISO/IEC 17025. The present list of methods reflects the amendments adopted by the 42nd Session of the Codex Alimentarius Commission in 2019.

<sup>1</sup> The most updated version of the method should be used, in application of ISO/IEC 17025. The present list of methods reflects the amendments adopted by the 42nd Session of the Codex Alimentarius Commission in 2019.

# Codex Standard 234-1999 Recommended Methods of Analysis and Sampling



- 2003 – publication by Enrique Méndez of R5 sandwich ELISA + cocktail extraction
- 2003 – international collaborative study with RIDASCREEN® Gliadin and Cocktail (patented)
- 2006 – Méndez method endorsed as Codex Alimentarius Type I Method
- 2008 – Méndez method became Type I Method and Codex Standard 118-1979 revision

## CSX 234-1999

### PART A – METHODS OF ANALYSIS BY COMMODITY CATEGORIES AND NAMES

Gluten-free foods	Gluten	Enzyme-Linked Immunoassay R5 Mendez (ELISA) Method <i>Eur J Gastroenterol Hepatol</i> 2003; 15: 465-474	Immunoassay	I
-------------------	--------	--	-------------	---

## 2.1.1 *Gluten-free foods*

Gluten-free foods are dietary foods

- a) consisting of or made only from one or more ingredients which do not contain wheat (i.e. all *Triticum* species, such as durum wheat, spelt, and khorasan wheat, which is also marketed under different trademarks such as KAMUT), rye, barley, oats<sup>1</sup> or their crossbred varieties, and the gluten level does not exceed 20 mg/kg in total, based on the food as sold or distributed to the consumer, and/or
- b) consisting of one or more ingredients from wheat (i.e. all *Triticum* species, such as durum wheat, spelt, and khorasan wheat, which is also marketed under different trademarks such as KAMUT), rye, barley, oats<sup>1</sup> or their crossbred varieties, which have been specially processed to remove gluten, and the gluten level does not exceed 20 mg/kg in total, based on the food as sold or distributed to the consumer.

## 5.2 Method for determination of gluten

Enzyme-linked Immunoassay (ELISA) R5 Mendez Method.

## Gluten Analysis in Codex Alimentarius

**Codex Alimentarius Type I method for gluten analysis in food**

Standard CXS 234-1999 (adoption from 2019)

**R5 ELISA**


RIDASCREEN® Gliadin



**Méndez cocktail**


for gluten extraction from  
heat-processed food samples  
→ Cocktail (patented)

# Most comprehensive gluten testing portfolio ...



## RIDA®QUICK – Lateral flow tests

- **RIDA®QUICK Gluten quant.**  
AOAC OMA approval in preparation
- RIDA®QUICK Gliadin  
AOAC-OMA 2015.16  
AOAC-PTM 101702



## RIDASCREEN® – ELISA

- RIDASCREEN® Gliadin  
**AOAC-OMA 2012.01**  
AOAC-PTM 120601
- **RIDASCREEN®EASY Gluten**  
AOAC PTM approval in preparation
- RIDASCREEN® Total Gluten  
AOAC-OMA 2018.15
- RIDASCREEN® Gliadin competitive  
AOAC-OMA 2015.05
- RIDASCREEN®FAST Gliadin sens.



## SureFood® – Real-time PCR

- SureFood® ALLERGEN Gluten
- SureFood® ALLERGEN 4plex Cereals / Getreide

## From farm to fork – gluten testing along the food production chain







## Grain harvest and processing



- Procurement and processing of raw materials
- Detection of contamination



RIDA®QUICK **New**  
Gluten quant.




RIDA®SMART BOX  
Alternative: smart phone  
(Android)



Data management

# Integrated gluten management




## The functioning of RIDA® SMART BOX

- 1 An image of test bands of lateral flow strip is taken
- 2 The data is transferred to the RIDA® SMART APP
- 3 The RIDA® SMART APP calculates the results
- 4 Full connectivity e.g. cloud, e-mail, pdf, excel, printer


# System solution

Customer A



Data export

Global QC management



Site C  
Data export


Site A  
Data export




Customer B



Data export



Site B  
Data export



Customer C



Dashboard with real-time data

Date & Time	Appl	Application	Batch Number	Customer ID	Notes	Product Code	Product Name	Result	Sample ID		
2023-12-06 15:42:30	allerge	ns	N/A	test		R7103	RIDA-QUICK Soya	N/A	5	null	SAS Systemintegrati on
2023-12-06 14:49:38	allerge	ns	Food Test	123455		RAL7073	RIDA-QUICK Gluten quant	< 2 mg/ kg	4	null	SAS Systemintegrati on
2023-12-06 12:48:58	allerge	ns	Food Test	123455		RAL7073	RIDA-QUICK Gluten quant	< 2 mg/ kg	3	null	SAS Systemintegrati on
2023-12-06 12:48:47	allerge	ns	Food Test	123455		RAL7073	RIDA-QUICK Gluten quant	< 40 mg/kg	2	null	SAS Systemintegrati on
2023-12-06 12:48:19	allerge	ns	Food Test	123455		RAL7073	RIDA-QUICK Gluten quant	< 2 mg/ kg	1	null	SAS Systemintegrati on
2023-12-06 15:35:32	allerge	ns	N/A	test	test	R7003	RIDA-QUICK Gladin	POSITIVE	test	null	SAS Systemintegrati on
2023-12-06 15:31:56	allerge	ns	N/A	test	1	R7003	RIDA-QUICK Gladin	N/A	1	null	SAS Systemintegrati on

- Dashboard with real-time data and analytics
- Customized APP design
- Data export to your own system via API





## Intermediate trade

Management of raw materials and ingredients for food industry



RIDA®QUICK **New**  
Gluten quant.



RIDA®SMART BOX

+



Data management

or



Samples to laboratory



## Quality management

Ingredients' handling  
Cross-contamination



RIDA®QUICK **New**  
Gluten quant.



RIDA®SMART BOX

+



Data management

or



RIDASCREEN® **New**  
EASY Gluten

+



Automation/  
LIMS



# Laboratory

Efficient sample handling



Fast, improved R5-ELISA method



**RIDASCREEN® EASY  
Gluten**

+



New easy extraction  
method included



ELISA automation

or



Other available ELISA  
test kits depending on  
customer needs



## Inline control production



- Cross-contamination during production
- Cleaning control



RIDA®QUICK **New**  
Gluten quant.



RIDA®SMART BOX

+



Data management



ready-to-swab

RIDA®QUICK Gliadin

or



Samples to laboratory



# Official control of foodstuffs

Consumer safety



RIDASCREEN®  
Gliadin

+



Automation/  
LIMS

- R5 ELISA with Cocktail (patented) – Méndez ELISA
- Codex Alimentarius type I reference method
- AOAC Official Method of Analysis for foods
- AOAC-PTM 101702

Confirmation/  
investigation




SureFood® ALLERGEN  
Gluten




## Most comprehensive gluten testing portfolio ...

... based on R5 antibody to ensure result comparability between different methods



### RIDA®QUICK – Lateral flow tests

- RIDA®QUICK Gluten quant.  
AOAC OMA approval in preparation
- RIDA®QUICK Gliadin  
AOAC-OMA 2015.16  
AOAC-PTM 101702



### RIDASCREEN® – ELISA

- RIDASCREEN® Gliadin  
AOAC-OMA 2012.01  
AOAC-PTM 120601
- RIDASCREEN®EASY Gluten  
AOAC PTM approval in preparation
- RIDASCREEN® Total Gluten  
AOAC-OMA 2018.15
- RIDASCREEN® Gliadin competitive  
AOAC-OMA 2015.05
- RIDASCREEN® Gliadin sensitive



### SureFood® – Real-time PCR

- SureFood® ALLERGEN Gluten

## Traceability due to same immunological target (R5)

r-biopharm®



18

Codex Type I Method  
(Méndez method)

RIDASCREEN® Gliadin  
+ Cocktail (patented)

RIDASCREEN® Gliadin  
+ Cocktail ECO

➤ more convenient and  
environment-friendly

RIDASCREEN® EASY Gluten  
+ Extraction tablet

➤ fast ELISA method and  
easy extraction procedure

RIDA® QUICK Gluten quant.  
+ Extraction tablet

➤ lateral flow method and  
easy extraction procedure

RIDASCREEN® Total Gluten  
+ Cocktail (patented)

➤ Specialty method for oats  
according to AOAC SMPR

# Quality management at laboratories

## Quality assurance

We recommend strongly


- to run standards and samples in duplicate
- to run control samples with every run  
(e.g. R7012 Set of 3 Processed Gliadin Assay Controls)

## Standard curve compliance criteria

- OD of Standard 6 > 1.2
- Lot-specific standard curve from CoA
- It is difficult to indicate exactly further criteria, since
  - each lab uses different equipment
  - each lab has its own environmental conditions (temperature, humidity, ...)
  - each lab has other people running the test

## How to use the CoA's standard curve information?

CoA	OD	CV	B/Bmax	Diff. between standards
Std. 1	0.131	3.1%	5.3%	
Std. 2	0.591	1.0%	24.1%	18.7%
Std. 3	0.970	3.1%	39.5%	15.4%
Std. 4	1.433	5.0%	58.4%	18.9%
Std. 5	2.027	5.1%	82.6%	24.2%
Std. 6	2.455	4.0%	100.0%	17.4%



- Your actual standard curve should have a similar shape; only one inflection point
- OD values can be different; B/Bmax and difference between standards should be similar
- CV of duplicates up to 15% are ok, but can be also higher for Standard 1

## We recommend

- to calculate mean value, SD, and CV from the first 10 runs of standards in duplicate (20 single OD values)
- to determine B/B<sub>max</sub> for standards and distance between standards from these runs as target and define a target range
- to run control samples (one, better two) with each run
- calculate mean values, SD, and CV as target from the first 10 runs and define a target range; plot values in a diagram

# Thank you for your attention!

More about R-Biopharm  
Food & Feed Analysis



<https://r-b.io/food>

