


PFAS Analysis | Building a Comprehensive Capability

Khalid Ghaffar
Regional Business Development Manager



Perfluoro and Polyfluoro alkylated Substances (PFAS)



- PFAS = Poly- or Per- fluorinated hydrocarbon chains with various functional groups on at least one end.
- First created in the 1930s with widespread applications
- Stable and persistent in the environment (POP) and Bio-accumulative
- Identified in environmental samples worldwide
- Most humans have PFAS in their blood





Dangerous PFAS Chemicals Are in Your Food Packaging

CR found 'forever chemicals' in bowls, bags, plates, and wrappers, even from some companies that say they've phased them out

PFAS contamination may be an emerging food safety issue
Food Business News (September 14, 2021)

C&EN “Pesticide breakdown products found in hundreds of US streams
Extensive environmental survey of pesticides and their transformation products reveals potential for hidden toxicity”



Sludge explained: What you need to know about “forever chemical” contamination on some Maine farms
By Annie Ropeik Maine
PUBLISHED 6:00 AM ET Mar. 08, 2022



Some of the worst PFAS contamination in Maine is linked to sludge from paper mills that used the chemicals in protective coatings, but almost all biosolids tend to contain at least low levels of PFAS.

CBS NEWS
**HEALTH
WATCH**

Raincoats, undies, school uniforms: Are your clothes dripping in "forever chemicals"?

healthwatch

BY HANNAH NORMAN

MARCH 30, 2023 /5:00 AM /KAISER HEALTH NEWS



HAWAII
NEWS NOW



Bottled water distributed to hundreds after toxic 'forever chemicals' found in Kunia well

Published: Mar. 24, 2023

The Guardian



Toxic 'forever chemicals' found in toilet paper around the World

Research finds waste flushed down toilets and sent to sewage plants probably responsible for significant source of water pollution

CBS NEWS

Eating one fish from U.S. lakes or rivers likened to drinking month's worth of contaminated water

JANUARY 17, 2023 / 8:04 AM / AFP



Why do we care?

Exposure to PFAS can lead to adverse health outcomes in humans



Increased cholesterol levels



Increased risk of high blood pressure or pre-eclampsia in pregnant women



Decreased vaccine response in children



Small decreases in infant birth weights



Changes in liver enzymes




Increased risk of kidney or testicular cancer



cancer (for PFOA), and thyroid hormone disruption (for PFOS, PFBS).

Harvard University, School of Public Health
“PFAS exposure linked with worse COVID-19 outcomes”

People can be exposed through water, air, soil, packaging, food, and more




Testing for PFAS requires a comprehensive analysis solution



Sample prep

- SPE cartridges
- Sample preparation apparatus (manual vacuum manifold or an automated system)
- Nitrogen sample evaporator




Instruments

- Ultra-Performance Liquid Chromatography (UPLC)
- Tandem quadrupole mass spectrometer (MS/MS)
- PFAS hardware kit
- Analytical columns



Informatics

- Instrument control and acquisition software
- Data processing and reporting
- Output to LIMS system or similar database
- Local quality system



CRM and PT

- Certified reference materials (CRM)
- Proficiency testing (PT) program for analysts

PFAS Testing Options – Total Solution for Successful Applications

Waters™

Consumables + Instrument + Methods + Software + Customer Education + Services

Xevo™
TQ-S micro



SPE enrichment prior to analysis.
Sample prep allows for use of mid-level sensitivity for MS

Xevo™
TQ-XS



Large volume/direct inject applications and SPE enrichment methods. Limited sample volume situations (e.g., serum). Applications requiring trace level detection limits

Xevo™
TQ ABSOLUTE



Achieve ultra-lower limits of quantitation to meet regulatory requirements with the most sensitive LC-MS/MS on the market

Chemistry

PFAS Kit with isolator column, analytical column, PEEK tubing, cartridges, standards, and more



Oasis WAX for PFAS is a quality controlled SPE



Sep-Pak PS2
6 cc Cartridge, for Method 537.1




Data Management





Getting Ready for Analysis: LC modifications for PFAS

No LC Modifications

PFOA standard




PFOA in blank




PFAS Kit Installed

PFOA standard




Background PFOA contaminant

No PFOA peak in blank



What is a Tandem Quadrupole MS?

- Tandem Quadrupoles consist of 2 mass analysers(quads) and a collision cell, enabling **more selective** analyses to be performed compared to a single quadrupole or other non MS detection mode.
- Their selectivity has made them the industry standard for MS based quantitation.



LC-MS/MS has been selected by regulators for PFAS for the technology's selectivity and sensitivity

Waters™

- High selectivity ➤ Reduce or eliminate matrix interferences
- High sensitivity ➤ Trace levels - low reporting limits
- Quantitative accuracy ➤ Reproducibility, stability, dynamic range and accurate quantitation of targets at low levels in matrix
- Robustness ➤ Complex sample matrices, reduced sample clean-up




Waters ACQUITY™ Premier UPLC™ + Xevo™ TQ Absolute MS



Waters ACQUITY UPLC + Xevo TQ-S micro MS

Multiple Reaction Monitoring (MRM)



Targeted PFAS Environmental Methods Using LC-MS/MS

Waters™

	<u>ASTM 7979 / 8421</u>	<u>ASTM 7968 / 8535</u>	<u>EPA 533</u>	<u>EPA 537.1</u>	<u>EPA 1633</u>	<u>EPA 8327</u>	<u>EU 2020/2184</u>	<u>EU 2022/2388</u>	<u>ISO 25101</u>	<u>ISO 21675</u>
Matrix	SW GW WW	Soils	DW	DW	Aqueous (not DW), Soil, Solids, Fish Tissue	Non potable water (SW, GW, WW)	DW	Food	DW GW SW	DW GW SW WW < 2 g/L solid particulate matter
Prep	Dilute, filter	Dilute, filter	Oasis WAX	Sep-Pak PS2	Oasis WAX GCB	Dilute, filter	Direct Inject	Oasis WAX	Oasis WAX	Oasis WAX
Waters Collateral	App Note 720006329	App Note 720006764	App Note 720006808	App Note 720006695	App Note 720008117	App Note 720006329	App Note 720007413	App Note 720007482	App Note 720006471	App Note 720006471
Column	CSH Phenyl Hexyl 2.1x100mm	CSH Phenyl Hexyl 2.1x100mm	BEH C18 2.1x100mm	BEH C18 2.1x100mm	BEH C18 2.1x50mm	CSH Phenyl Hexyl 2.1x100mm	ACQUITY Premier BEH Shield RP18 Column	BEH C18 2.1x100mm	BEH C18 2.1x100mm	BEH C18 2.1x100mm
Minimum Instrument	TQ Absolute	TQ Absolute	TQ-S micro	TQ-S micro	TQ-S micro	TQ Absolute	TQ Absolute	TQ Absolute	TQ-S micro	TQ-S micro
Number of compounds	24 / 44	24 / 44	25	18	40	24	20	30	2	30



On Apr 10, 2024 US EPA finalized **historic regulations** to limit PFAS in drinking water

PFAS National Primary Drinking Water Regulation overview:

- An enforceable 4 parts per trillion (ppt) limit on PFOA and PFOS in drinking water
- A non-enforceable maximum contaminant level goal for PFOA and PFOS at zero
- A limit of 10 ppt on three other categories of PFAS: PFNA, PFHxS, and “GenX” chemicals
- A hazard index for mixtures of 2+ of the above PFAS, plus PFBS
- Public Water Systems (PWS) must complete initial monitoring by 2027 and implement solutions by 2029. From 2027+, they will be required to complete ongoing compliance monitoring

Announcement of the new regulations here:
<https://www.epa.gov/sdwa/and-polyfluoroalkyl-substances-pfas>


Putting ng/L into perspective

Nanograms of compound per liter of solution (aka part-per-trillion (ppt))

– **One** part per 1,000,000,000,000 (10^{12}) parts

Some ways to visualize one ng/L


1 drop of water in 20 Olympic-sized swimming pools



One second of time in approximately 31,700 years



Traveling 6 inches out of a 93-million-mile journey to the sun¹



1. The earth is 93 million miles (avg.) from the sun

Analytical Requirements: Drinking Water Methods Required to Support Testing

Method 537.1: Determination of Selected Per- and Polyfluorinated Alkyl Substances in Drinking Water by Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

HFPO-DA	PFNA
NEtFOSAA	PFOS
NMeFOSAA	PFOA
PFBS	PFTA
PFDA	PFTTrDA
PFDoA	PFUnA
PFHpA	11Cl-PF3OUds
PFHxS	9Cl-PF3ONS
PFHxA	ADONA



SepPak PS2 cartridge

Method 533: Determination of Per- and Polyfluoroalkyl Substances in Drinking Water by Isotope Dilution Anion Exchange Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry

11Cl-PF3OUds	4:2FTS
9Cl-PF3ONS	PFHxS
ADONA	PFHxA
HFPO-DA	PFMPA
NFDHA	PFMBA
PFBA	PFNA
PFBS	6:2FTS
8:2FTS	PFOS
PFDA	PFOA
PFDoA	PFPeA
PFEESA	PFPeS
PFHpS	PFUnA
PFHpA	



Oasis WAX cartridge

537.1: Waters App Note [720006695](#)

533: Waters App Note [720006808](#)

Targeted Analysis with Enhanced Sensitivity



Xevo TQ Absolute

Xevo TQ Absolute - Key Sustainability Benefits

Harness Absolute Power



XEVO[™]
TQ ABSOLUTE

- Up to **15x more sensitivity** for challenging negative ionizing compounds due to detector design improvements. Xevo TQ Absolute **Signal:Noise >5 million to one** for reserpine and chloramphenicol compared to **Agilent's 6495D, 4 million to one**.
- Uses **45% less bench space** than the leading competitor
- Uses **~50% less power & gas** and generate **~50% less heat** than the leading competitor
- Save up to **50% on data review** time using `waters_connect` for Quantitation
- **MyGreenLab certification:** Xevo TQ Absolute scores well on a range of sustainability criteria including manufacturing, shipping, product and packaging content, end of life and energy consumption.
- **Environmental Impact Factor:** Xevo TQ Absolute consumes **one third of the in-use energy** compared to **Agilent 6495D** and the **overall ACT label* rating score was half that of the Agilent 6495D**.

*The **ACT Environmental Impact Factor Label** was designed to address the need of both scientists and procurement specialists for clear, third-party verified information about the environmental impact of laboratory products.

Waters[™]

		US
Waters Xevo TQ Absolute Triple Quadrupole Mass Spectrometer		
Wexford, Ireland SKU 176005117		
Environmental Impact Scale Decreasing Environmental Impact		
1 10		
Manufacturing		
Manufacturing Impact Reduction		1.0
Renewable Energy Use		Yes
Responsible Chemical Management		1.0
Shipping Impact		7.9
Product Content		5.0
Packaging Content		5.0
User Impact		
Energy Consumption (kWh/day)		19.3
Water Consumption (gallons/day)		N/A
Product Lifetime		3.0
End of Life		
Packaging		5.2
Product		6.5
Innovation		
Innovative Practices		-1.0
Environmental Impact Factor:		52.9
Label Valid Through:		December 2025
		act.mygreenlab.org

June 2023



EPA Announces New Drinking Water Health Advisories for PFAS Chemicals, \$1 Billion in Bipartisan Infrastructure Law Funding to Strengthen Health Protections

Compound	US EPA Health Advisory Level (ng/L)	UCMR 5 Minimum Reporting Level (ng/L)
PFOA	0.004 (interim)	4
PFOS	0.02 (interim)	4
PFBS	2000	3
GenX	10	5

- PFAS from fluoropolymers and coating are ubiquitous in common sampling and analytical equipment and can easily contaminate solvents and reagents
- Laboratory contamination is common and near impossible to eliminate
 - Perfluorobutane sulfonic acid and its potassium salt (PFBS) extremely common
- Such contamination leads to reporting of false positives and/or higher reporting limits due to elevated background

Essential mitigation options:

- Avoid all products likely to contain fluoropolymers – e.g. vials with Teflon seals. Use HDPE, polypropylene or polyethylene materials
- Minimize risks from contaminated dust and air
- Screen all analytical consumables, solvents and reagents including water purifying systems
- Carry out extensive use of procedural blanks from sample storage
- Traceability systems for the consumables used in the laboratory
- Training !
- Modifications to the LC system fitting a PFAS Installation Kit

Sensitivity – Controlling the Background in a Routine Laboratory


[WHITE PAPER]

Waters™

Best practices for monitoring PFAS contamination in a routine shared-space commercial laboratory

Nicola Dreolin, Henry Foddy, Kari Organtini, Stuart Adams, Ken Rosnack, Peter Hancock
Waters Corporation, Milford, MA, USA and Waters Corporation, Wilmslow, UK.

[Link to White Paper](#)




Common Sources of PFAS Contamination to Avoid

External Sources	Direct Sources
Clothing/Lab coats treated with waterproofing materials	PTFE (Teflon) containers lined caps, and tubing
Waterproof papers, notebooks, binders	Aluminum foil
Cosmetics and personal care products (sanitizers, lotions, etc.)	Pipette tips branded as being "low retention"
Teflon tape	Permanent markers
Latex gloves	Vacuum grease
Antifog eyewear wipes and sprays	Glass transfer pipettes
Soaps and dishwashing detergents	PTFE filters

What about SPE + Xevo TQ Absolute for EPA HALs?

Waters™




- Enrichment needed
A 500-fold sample enrichment was performed using Oasis™ WAX for PFAS SPE cartridges

250 ml to final volume 0.5 ml

- Highly Sensitive MS/MS
Enhanced negative ion mode sensitivity of the Xevo TQ Absolute
- Blanks will be challenging!
Background PFAS contamination from routine materials in laboratories must be highly controlled


[Best practices for monitoring PFAS contamination in a routine shared-space commercial laboratory](#)

Sensitivity - EPA HALs Testing



- It may be impossible to completely eliminate contamination during sample preparation but it can be controlled to an acceptable amount.

- Critical assessment of the background of the laboratory vs the response from the spiked samples
- Blanks response <30% of the method LOQ




Limit of Quantitation at 0.001 ng/l (ppt)

Compound	Labratory 1 LOQ ng/L	Labratory 1 S:N of LOQ	Labratory 2 LOQ ng/L	Labratory 2 S:N of LOQ	EPA HAL ng/L	Interim minimum reporting level ng/L
PFOA	0.001	10	0.001	10	0.004 (interim)	4
PFOS	0.001	12	0.001	27	0.02 (interim)	4
PFBS	0.0009	24	0.002	27	2,000	3
HFPO-DA (GenX)	0.004	21	0.004	10	10	5

Limit of quantitation (LOQ) and signal:noise (S:N) for each PFAS determined in both laboratories using two different TQ Absolute systems

EPA method 1633 Available for Non-Drinking Water Matrices:



Newly finalized method (January 2024) to test for **40 PFAS compounds** in **wastewater, surface water, groundwater, soil, biosolids, sediment, landfill leachate, and fish tissue**

Created by the **US EPA's Office of Water**, partnership with the **US Department of Defense's Strategic Environmental Research and Development Program (SERDP)**

Validated through a **multi-laboratory study** to finalize method and add performance criteria

Intended to be used for:



Clean Water Act Compliance (wastewater discharge permits – NPDES)




Superfund sites



Remediation and investigation programs (including for DoD installations)

Demonstration of Performance for 1633




46 native PFAS covering both EPA 1633
and ASTM 8421 lists


- ❑ Minimum Levels on the TQ Absolute are 20x lower than those reported in the EPA 1633 draft methods (on the Xevo TQ-S micro)
- ❑ S:N values at the minimum level were all ≥ 3 as required by EPA 1633

All initial calibration requirements in EPA 1633 are easily met on the Xevo TQ Absolute

LC Gradient Optimization




- Cholic acids (TDCA, TCDCA, TUDCA) can cause interference with PFOS in the mass spectrometer due to the similarity in parent and fragment masses
- EPA 1633 requires a 1 min resolution of the cholic acids from PFOS
- Methanol is typically used for PFAS methods, but does not resolve cholic acids from PFOS
- Acetonitrile provides more than 1 minute separation of cholic acids from all PFOS isomers



PFAS analysis by EPA 1633 – results from LC-MS/MS


Certified reference material (CRM) validation

PFAS in wastewater CRM



Mean trueness of 92%
Trueness range of 73 – 112%




PFAS in soil CRM



Mean trueness of 97%
Trueness range of 85 – 120%

Solid-phase extraction (SPE) - why we use it for PFAS analysis

- **Regulatory compliance**
 - US EPA 533, 537.1, and 1633 require solid phase extraction clean up steps
 - Reproducibility across laboratories running same method
- **Clean up from complex matrices**
 - Environmental waters, solids, biosolids, and other samples contain interferences
 - Reduce ion suppression and background
- **Concentration of PFAS**
 - Measuring at low levels (ng/L)
 - Increase limits of detection
- **Selectivity in extraction**
 - Use of WAX sorbent more specific to PFAS



The Guardian

Thousands of toxins from food packaging found in humans - research

Metals and PFAS linked to serious health issues are among compounds found, highlighting need for further scrutiny



Research found that among the worst offenders is plastic, a material that is largely unregulated and can contain thousands of chemicals. Photograph: Aleaimage/Getty Images/iStockphoto

Tom Perkins
Fri 27 Sep 2024 07.00 EDT

'Forever Chemicals' Are Found in Some Milk, Including Organic


A Consumer Reports investigation highlights gaps in how the U.S. tests and regulates PFAS in food



PHOTO ILLUSTRATION: CHRIS GRIGGS/CONSUMER REPORTS, GETTY IMAGES

By Lauren Kirchner · Data visualizations by [Andy Bergmann](#)
May 2, 2024

Sample Preparation Protocols for Food



Produce and Milk





Produce	Commodity Class
Romaine Lettuce	High Water
Russet Potato	High Starch
Whole Carrot	Low Water
Strawberry	High Sugar
Cranberry	High acid and sugar

Produce: Waters App Note [720007333](#)

Milk App Note coming soon

QuEChERS Method for Produce and Milk


Waters™



**Quick, Easy,
Cheap, Effective,
Rugged, Safe =
QuEChERS**

Extracted Ion Chromatograms


Waters™




ACQUITY UPLC™ I-Class
with Xevo™ TQ-XS

Linear and Branched Isomers for Complete Assessment


Cranberry



Potato



Carrot



PFAS Quantitation in Milk



1.

5 g produce
Spike Surrogates at 1 ng/g (MPFAC-24ES + M3HFPO-DA)
Add 5 mL of water

2.

Add 10 mL Acetonitrile + 150 μ L Formic Acid
Shake for 1 minute

3.

Add AOAC QuEChERS Salts (6 g $MgSO_4$, 1.5 g Sodium Acetate)
Shake 5 minutes
Centrifuge 5 minutes at 4,000 rpm

4.


Transfer 5 mL supernatant to 15 mL dSPE tube (900 mg $MgSO_4$, 300 mg PSA)
Shake 5 minutes
Centrifuge 5 minutes at 4,000 rpm

5.

Dilute 1:1 in 2mM ammonium acetate
Centrifuge 5 minutes at 4,000 rpm
Spike Internal Standard equivalent to 1 ng/g (MPFAC-C-IS)

Milk App Note [720007687](#)

Extracted Ion Chromatograms 0.1 ng/g each PFAS




Matrix Matching or Isotope Dilution?

Matrix Matched Curve

- Use only two internal standards
- $^{13}\text{C}_8$ -PFOA and $^{13}\text{C}_8$ -PFOS
- Makes an estimate for recovery
- Calibration is performed in similar conditions as sample

Solvent Curve


- Use isotope dilution
- 23 ^{13}C and D_2 labelled standards used
- Accounts more accurately for recovery
- Internal standards also account for matrix effects




Peak area comparison of PFHxS, PFTreDA, GenX, and NMeFOSAA in 100 ng/L points of both the solvent and matrix matched curves.

Recovery of PFAS in Spiked Milk Samples - What Works Best?

Recovery of 0.10 ng/g spike




Recovery of 1.0 ng/g spike



Mean percent accuracy





- ❖ 85% matrix matched - less costly (stds) but need an available clean matrix; how many different matrices being supported?
- ❖ 97% solvent curve with isotope dilution – can be expensive to obtain all the standards, however, greater accuracy and more flexible for wider range of sample types

Sample Type
Whole Egg
Salmon
Tilapia
Ground Beef
Beef Liver
Beef Kidney



Waters App Note [720007482](#)

Salmon Spiked at 0.1 µg/kg – Digestion & WAX SPE Prep




ACQUITY BEH C₁₈ (2.1 x 100 mm, 1.7 µm)

- Chromatograms shows the 30 PFAS evaluated in the method spiked into salmon extract at 0.1 µg/kg (100 ppt)
 - Peak assignments are located in the Appendix of the [application note](#)


Even though the extraction method results in 5x dilution of sample concentration, the sensitivity of the TQ-XS allows for trace level detection of PFAS in matrix samples.

Taurodeoxycholic acid (TDCA) interference in meat

(A) chromatograms of PFOS MRM transitions 499 > 91 and 499 > 80 and M8-PFOS in beef liver, chicken egg, and ground beef showing taurodeoxycholic acid (TDCA) as an interferent in 499>80 transition in liver and egg.



(B) M-PFOS injection standard peak area (y-axis) over 12 sample injections of each sample matrix (x-axis): salmon, tilapia, beef, liver, and kidney



Cosmetics

Waters™



1 g sample
+ 10 mL
methanol

Sonicate 30
minutes

Filter with
0.22 um
GHP syringe
filter

Dilute 1:1
with 2mM
ammonium
acetate

Sample	compound	conc [ng/mL]	conc [ng/g]
Eye Shadow	PFBA	0.29	2.9
Eye Shadow	PFHxA	0.02	0.2
Eye Shadow	PFOA	0.12	1.2
Eye Shadow	PFTeDA	0.02	0.2
Eye Shadow	PFBS	0.02	0.2
Foundation	PFBA	0.18	1.8
Foundation	PFOA	0.01	0.1
Foundation	PFTeDA	0.02	0.2

[Environ. Sci.: Processes Impacts, 2018, 20, 1680](#)

96 Well WAX SPE procedure for Human Serum and Plasma

Waters™

Pretreat Samples (50 µL) – spike isotope dilution standards, PPT with 3x ACN, dilute with 5x 1% Formic Acid

Condition Plate

200 µL 2% ammonia/methanol
200 µL methanol
200 µL 1% Formic Acid

Load Sample

Wash Plate

200 µL 1% Formic Acid
200 µL 1:3 ACN:Water

Elute

50 µL methanol
50 µL 2% (v/v) ammonia/methanol → 2x
Dilute with 50 µL 2% ammonium acetate containing Injection Standards




Oasis WAX 96 well
µElution plate

2 mg per well

Waters App Note [720007114](#)


Human Pooled Serum and Plasma Samples

Serum Lot 1




Total PFAS =
18.96 ng/L

Serum Lot 2




Total PFAS =
14.44 ng/L

Serum Lot 3



Total PFAS =
8.86 ng/L

Plasma



Total PFAS =
1.71 ng/L

Education and training are fundamental to successfully establishing a PFAS analysis capability

Waters can help!

- Custom onsite training
- Virtual training
- Advanced method consultation
- Education subscriptions
- Certified in-house courses
- E-learning courses



No matter what you choose, we support you.

The Waters™ Customer Education Team is here for you with expert instructors whether you have just purchased a new Waters instrument, need help learning a new software or application, or simply require refresher training on any of our platforms.

A Comprehensive Workflow Solution for PFAS

Waters™



Oasis SPE




PFAS Base Hardware Only Kit



Atlantis Premier BEH C18 AX 50 mm used as Isolator Column




waters_connect for Quantitation



ACQUITY Premier Analytical Columns
Atlantis Premier AX Analytical Column

Xevo TQ Absolute

ACQUITY Premier



Analytical Professional Services



Native PFAS Precision and Recovery Standard Solution



PP Vials & Caps

Waters™ | ERA

Certified Reference Materials

Quick Response PTs

Proficiency Tests




Waters™

Waters



Beyond PFAS – LC-MS/MS Applications for Environmental Testing Waters™

PPCPs in Water using LC-MS/MS



SOLUTIONS

- Targeted method for the determination of PPCPs in wastewater using US EPA Method 1694.
- UPLC with Tandem Mass Spectrometry water
- Suitable for monitoring various water sample matrices with minimal sample preparation
- 200 mg Oasis HLB Solid Phase Extraction (SPE) Cartridges


Basic UPLC conditions at 100 ng/L, fortified in surface water

Peak	UPLC	Class	RT	MS/MS (m/z)	Waters	MS/MS (m/z)
Chlorhexidine	MacroSIL	0.001	1.1	50	50	21
Triclosan	MacroSIL	0.010	0.5	50	110	21
Propylparaben	MacroSIL	1.000	3.00	61	78	—
Propylparaben d13	MacroSIL	1.000	3.00	61	110	—
Erythromycin	MacroSIL	1.000	3.00	50	50	21
Amoxicillin	MacroSIL	0.700	1.1	50	50	21
Hydrocortisone	Tandem	0.001	1.1	75	100	—
Cephalosporin	Phenomenex	0.000	0.20	50	75	—
Sulfamethoxazole	Phenomenex	0.000	0.20	60	62	—
Indinavir sulfate	Phenomenex	0.000	0.1	50	50	21

These macrolides and six other important antibiotics analyzed with off-line SPE and UPLC-ESI-MS/MS in positive mode at basic and acidic UPLC conditions

8001 Waters Corporation COMPANY CONFIDENTIAL

Chlorinated Acid Herbicides in Water



Analyze, Without the Need to Derivatize

CHALLENGES


- Lengthy preparation steps with derivatization
- Poor chromatography with historic GC-ECD, or GC-MS methods

SOLUTIONS

- Targeted method for the determination of 20 acidic herbicides in water samples that is suitable for monitoring both drinking water and surface waters with minimal sample preparation
- No lengthy and dangerous derivatization required
- Short analytical run time

8001 Waters Corporation COMPANY CONFIDENTIAL

Link to Temple University Publication on PPCPs in the Delaware River



Abstract

Occurrence and aquatic toxicity of contaminants of emerging concern (CECs) in tributaries of an urbanized section of the Delaware River

Waterford

Keywords: Waters™, Capillary Electrophoresis, Robert Harnack, Jonathan Hart, A. Ronald McElroy


Abstract: The presence of contaminants of emerging concern (CECs) in environmental matrices is an ongoing issue. The present project was carried out to increase our understanding of the loading, distribution and potential risk (CEC) by sampling large and small tributaries in a specific area of the Delaware River watershed to surface water. Data for a highly chemical and significantly impacted waterway measure plan (WMP). Fifteen target compounds were selected for analysis based on their high frequency of detection in a previous multi-media study conducted on the Delaware River estuary. Ten sampling sites were chosen on tributaries receiving untreated municipal and industrial discharges. Sampling locations were above and below potential source discharges. Sampling was designed to assess seasonal differences in CECs loading. The seasonal environmental concentrations of the target compound present a detailed view of their spatial and temporal impacts on urbanized receiving water. An index of source loading system was applied to the sample location to compare temporal and seasonal concentrations, using target compound ratio profiles relative to predicted no effect levels and detecting a common seasonal variable. Toxicological and bioaccumulation assessment of the compounds of high concern and CECs in the aquatic life of the Delaware River tributaries.

Table 1. Compounds of interest

Active	Therapeutic category	Description
Carbamazepine	Anti-epileptic drugs	Prescription drug, anticonvulsant, mood-stabilizer, treatment of certain deficit hyperactivity disorder (ADHD)
Clarithromycin	Antibiotic	Prescription drug, a macrolide antibiotic
Dihydrochloride	A metabolite of a nitroglycerin	An inactive metabolite of prescription cardiac drug nitroglycerin
Diltiazem	Anti-hypertensive	Prescription drug, anticonvulsant, mood-stabilizer, treatment of ADHD
Dihydrochloride	Antibiotic	Prescription and over-the-counter drug, an antibiotic used to treat otitis
Erythromycin	Antibiotic	Prescription drug, a macrolide antibiotic
Glimepiride	Anti-diabetic	Prescription drug, lipid-lowering agent
Guanfacine	Drug of abuse	Mefenamic acid
Ibuprofen	Non-steroidal anti-inflammatory	Prescription and over-the-counter drug, a non-steroidal anti-inflammatory drug used for pain
Metformin	Anti-diabetic	A prescription drug for diabetes treatment
Ranitidine	Antihistamine	Prescription and over-the-counter drug for peptic ulcer and gastroesophageal reflux disease treatment
Sulfamethoxazole	Antibiotic	Prescription drug, a bacteriostatic antibiotic frequently used in combination with trimethoprim
Thiamethoxam	Fungicide; pesticide	Prescription and veterinary drug in a parasiticide also used as an agricultural fungicide
Triclosan	Personal care product	Antibacterial, antiseptic agent
Trimethoprim	Antibacterial	Prescription drug, a bacteriostatic antibiotic frequently used in combination with sulfamethoxazole

8001 Waters Corporation COMPANY CONFIDENTIAL

Anionic Polar Pesticide Solution



Analyze, Without the Need to Derivatize

CHALLENGES

- Difficulty with chromatographic retention
- Several polar anionic pesticides share mass transitions or can degrade in solution or via in source fragmentation, generating interferences in MRMs and thus false positives.

SOLUTIONS

- Anionic Polar Pesticide Column exhibits stable retention times, even across different commodities.
- AMPA is the first eluting peak with retention 3.5 times greater than the column void volume.
- Column resolves such false detections

8001 Waters Corporation COMPANY CONFIDENTIAL

Want to learn more? Additional resources for PFAS in Food

- [Analysis of 28 EU Regulated and Recommended PFAS in Food via LC-MS/MS – Part 1: Vegetable, Fruit, and Baby Food](#)
- [Analysis of Per- and Polyfluoroalkyl Substances \(PFAS\) in Animal Products with an Enhanced Sensitivity LC-MS/MS Method using Fish Reference Materials as a Case Study](#)
- [QuEChERS Extraction of Per- and Polyfluoroalkyl Substances \(PFAS\) from Edible Produce with Sensitive Analysis on Xevo TQ-XS](#)
- [Total Workflow for the Sensitive Analysis of Per- and Polyfluoroalkyl Substances \(PFAS\) in Fish, Meat, Edible Offal, and Eggs](#)
- [Evaluation of the Performance of a Total Workflow Approach for the Analysis of Per- and Polyfluoroalkyl Substances \(PFAS\) in Fish, Using an Interlaboratory Study](#)
- [Matrix Matching or Isotope Dilution? A Comparison of Two Quantitation Approaches to Determine PFAS in Dairy Milk](#)
- [Optimization of Source and Transmission Parameters for a Mix of Labile and Stable Per – Or Polyfluoroalkyl Substances \(PFAS\) Using the Xevo™ G3 QToF Mass Spectrometer](#)

References

Company Name:	Eurofins Eaton Lab
Contact Person:	Dr. Bruce (Yongtao) Li
Telephone Number:	574-233-4777
Email:	yongtaoli@eurofinsus.com
Address Line 1:	110 S Hill Street
Address Line 2:	South Bend, Indiana 46617, USA
Comments	Lab has TQ-S, TQ-XS and TQ Absolute

Company Name:	Rhode Island State Health Laboratories
Contact Person:	Evan Philo, Principal Laboratory Scientist / Food Testing Coordinator
Telephone Number:	401-222-5553
Email:	Evan.Philos@health.ri.gov
Address Line 1:	50 Orms St.
Address Line 2:	Providence RI, 02904
Comments	TQ-Absolute for PFAS analysis in water

Company Name:	US EPA Region 5 Lab
Contact Person:	Larry Zintek
Telephone Number:	312 886 2995
Email:	Zintek.lawrence@epa.gov
Address Line 1:	77 West Jackson Blvd
Address Line 2:	Chicago, IL 60604
Comments	Multiple TQ-XS and TQ-Absolute systems for PFAS analysis

Company Name:	Suffolk County Water Authority
Contact Person:	Amanda Comando
Telephone Number:	631 218 1128
Email:	acomando@scwa.com
Address Line 1:	260 Motor Parkway
Address Line 2:	Hauppauge, NY 11788
Comments	Multiple TQ-S/TQ-XS systems

Company Name:	Complete Environmental Testing
Contact Person:	Jeff Smith, Lab Director
Telephone Number:	203 377 9984
Email:	jsmith@cetlabs.com
Address Line 1:	80 Lupes Dr
Address Line 2:	Stratford, CT 06615
Comments	Multiple TQ-XS for PFAS analysis

Company Name:	Katahdin Analytical Services
Contact Person:	Michael Flanders, President and CEO
Telephone Number:	207-874-2400
Email:	mflanders@katahdinlab.com
Address Line 1:	600 Technology Way
Address Line 2:	Scarborough, ME-04074
Comments	3 x TQ Absolute systems for PFAS analysis in drinking water and other matrices

Thank you

khalid_ghaffar@waters.com

