




# Military Laboratories For Quality Control MLQC




# INTRODUCTION

IN THE PURSUIT OF OPTIMIZING LIVESTOCK PRODUCTION, ZILPATEROL AND RACTOPAMINE HAVE EMERGED AS POWERFUL BETA-AGONIST DRUGS WIDELY USED AS GROWTH PROMOTERS IN CATTLE. THESE COMPOUNDS ENHANCE FEED EFFICIENCY, INCREASE LEAN MUSCLE MASS, AND REDUCE FAT DEPOSITION, MAKING THEM VALUABLE TOOLS IN THE MEAT PRODUCTION INDUSTRY. HOWEVER, THEIR USE IS NOT WITHOUT CONTROVERSY. WHILE APPROVED IN THE UNITED STATES AND CANADA, THESE DRUGS ARE STRICTLY BANNED IN THE EUROPEAN UNION AND MANY OTHER COUNTRIES DUE TO CONCERNS OVER POTENTIAL HEALTH RISKS TO CONSUMERS AND ANIMAL WELFARE ISSUES.

THE PRESENCE OF RESIDUES FROM THESE DRUGS IN ANIMAL TISSUES, SUCH AS MUSCLE, LIVER, AND KIDNEY, POSES A SIGNIFICANT CHALLENGE TO FOOD SAFETY. TO PROTECT PUBLIC HEALTH, IT IS ESSENTIAL TO MONITOR AND CONTROL THE LEVELS OF THESE RESIDUES, ENSURING THEY REMAIN WITHIN PERMISSIBLE REGULATORY LIMITS. THIS IS WHERE ADVANCED ANALYTICAL TECHNIQUES, SUCH AS LIQUID CHROMATOGRAPHY-TANDEM MASS SPECTROMETRY (LC-MS/MS), PLAY A CRITICAL ROLE. LC-MS/MS IS A HIGHLY SENSITIVE AND PRECISE METHOD FOR DETECTING AND QUANTIFYING TRACE AMOUNTS OF ZILPATEROL AND RACTOPAMINE IN COMPLEX MATRICES LIKE ANIMAL TISSUES.



# WHAT ARE ZILPATEROL AND RACTOPAMINE?





- **ZILPATEROL:**

- A BETA-AGONIST DRUG USED TO PROMOTE LEAN MUSCLE GROWTH IN CATTLE.
- APPROVED IN THE USA AND CANADA WITH SPECIFIC MAXIMUM RESIDUE LIMITS (MRLS).

- **RACTOPAMINE:**

- ANOTHER BETA-AGONIST USED TO IMPROVE FEED EFFICIENCY AND INCREASE LEAN MEAT PRODUCTION.
- ALSO APPROVED IN THE USA AND CANADA WITH DEFINED MRLS.

- **HEALTH CONCERNS:**

- POTENTIAL RISKS TO HUMAN HEALTH IF RESIDUES EXCEED PERMISSIBLE LIMITS.
- NEED FOR RELIABLE ANALYTICAL METHODS TO MONITOR RESIDUES IN FOOD PRODUCTS.



## PERMISSIBLE LIMITS IN THE USA, CANADA AND CODEX STANDARD:

### USA:

#### RACTOPAMINE:

MUSCLE: 30 MJG/G (PPB)

#### ZILPATEROL:

MUSCLE: 12 MJG/G (PPB)

### CANADA:

#### RACTOPAMINE:

MUSCLE: 40 MJG/G (PPB)

#### ZILPATEROL:

MUSCLE: 5 MJG/G (PPB)

### MAXIMUM PERMISSIBLE LIMIT IN CODEX:

#### RACTOPAMINE:

MUSCLE: 10 MJG/G (PPB)

#### ZILPATEROL:

MUSCLE: 0.5 MJG/G (PPB)

CODEX  
ALIMENTARIUS




Food and Agriculture  
Organization of the  
United Nations





World Health  
Organization

# ANALYSIS OF RACTOPAMINE AND ZILPATEROL BY USING UPLC-MS/MS





WE DETECT AND QUANTIFY RACTOPAMINE AND ZILPATEROL IN KIDNEY, MUSCLE, AND LIVER TISSUES TO ENSURE COMPLIANCE WITH REGULATORY LIMITS, USING UPLC-MS/MS AND ACCORDING TO OASIS PRIME MCX METHOD .



## Sample Collection and Preparation:

- \* Collect **5 g** of the sample (from liver, muscle, or kidney tissues).
- \* Ensure proper **homogenization** of the sample for consistency.

**Extraction:**  
Perform **extraction** to isolate the target analytes (ractopamine and zilpaterol) from the sample matrix by methanol.

**Solid-Phase Extraction (SPE):**  
Purify the extracted sample using the Oasis PRIME MCX SPE method to remove impurities and concentrate the analytes

**Sample Injection:**  
Inject the purified sample into the UPLC-MS/MS system for analysis.

**Detection and Quantification**





# WHAT WE HAVE DONE AND WHAT WILL WE DO?



**Analysis of Ractopamine and Zilpaterol in Cattle Tissues: From Method Development to Routine Application in Military Laboratories For Quality Control (MLQC)**




METHOD  
DEVELOPMENT



METHOD VALIDATION



IMPLEMENTATION OF  
THE METHOD



VERIFICATION  
AND ROUTINE  
ANALYSIS



# WHAT WE HAVE ACHIEVED IN MILITARY LABORATORIES FOR QUALITY CONTROL AND WHAT WAS THE RESULT ?








# RESULTS OF METHOD DEVELOPMENT THAT APPLIED IN MLQC

## LINEARITY:

- EXCELLENT LINEARITY WAS OBSERVED FOR BOTH RACTOPAMINE AND ZILPATEROL.
- CORRELATION COEFFICIENT ( $R^2$ ):**  $> 0.999$  FOR BOTH COMPOUNDS.





- **SELECTIVITY:**

- HIGH SELECTIVITY WAS ACHIEVED WITH WELL-RESOLVED PEAKS FOR RACTOPAMINE AND ZILPATEROL.
- NO INTERFERENCE FROM MATRIX COMPONENTS WAS OBSERVED.

- **CHROMATOGRAPHIC PEAKS:**



- **RACTOPAMINE:** RETENTION TIME = 3.15 MIN.
- **ZILPATEROL:** RETENTION TIME = 1.8 MIN.
- PEAKS WERE SHARP AND SYMMETRICAL, INDICATING GOOD SEPARATION.

Targetlyn XS - unritted

File Edit View Display Processing Window Help

### Ractopamine

Name	ID	RT	Area	Std. Conc	ppb	%Rec	S/N
R_Z_6_29_1_2025_01	6	3.10	401890	24.000	23.535	98.1	15267
R_Z_5_29_1_2025_01	5	3.10	308038	18.000	18.001	100.0	11493
R_Z_4_29_1_2025_01	4	3.10	245463	14.000	14.311	102.2	8290
R_Z_3_29_1_2025_01	3	3.10	141148	8.000	8.180	102.0	5233
R_Z_2_29_1_2025_01	2	3.10	37185	2.000	2.030	101.5	1382
R_Z_1_29_1_2025_01	1	3.09	19005	1.000	0.962	96.2	667
BLANK 29_1_2025_02	0	3.10	-19				1









### Zilpaterol (1)

Name	ID	RT	Area	Std. Conc.	ppb	%Rec.	S/N
R_Z 6_29_1_2025_01	6	1.77	681	0.960	0.886	98.5	70
R_Z 5_29_1_2025_01	5	1.77	545	0.700	0.714	102.0	43
R_Z 4_29_1_2025_01	4	1.76	378	0.500	0.500	100.0	38
R_Z 3_29_1_2025_01	3	1.77	220	0.300	0.299	99.8	18
R_Z 2_29_1_2025_01	2	1.73	65	0.100	0.102	102.4	5
R_Z 1_29_1_2025_01	1	1.75	16	0.040	0.038	97.6	2
BLANK 29_1_2025_02	0	1.94	8		0.030		1





# OUR NEXT STEP WITH UPLC-MS/MS

- WE WILL PERFORM **METHOD VALIDATION** (THE NEXT STEP) BY APPLYING THE FOLLOWING VALIDATION PROTOCOL

## •Preparation of Calibration Standards:

Prepare a series of standards with known concentrations of Ractopamine and Zilpaterol.

## •Spiked Samples:

Spike blank tissue samples with known concentrations of analytes for recovery studies

## Analysis by UPLC-MS

•Data Analysis:  
Calculate LOD, LOQ, precision, accuracy, and other parameters.



# ASPIRATIONS

The image features a dark blue background with several realistic water droplets of various sizes scattered around the edges. The droplets are rendered with highlights and shadows, giving them a three-dimensional appearance. The word "ASPIRATIONS" is centered in a bold, yellow, sans-serif font.

- **TO BECOME A GLOBALLY RECOGNIZED LEADER IN THE ANALYSIS OF VETERINARY DRUG RESIDUES AND FOOD ADDITIVES.**
- **ENSURING FOOD QUALITY AND SAFETY BY PROVIDING ACCURATE AND RELIABLE ANALYTICAL SERVICES.**
- **RECEIVING SAMPLES.**
- **OBTAINING ISO 17025 ACCREDITATION**

**THANK YOU FOR LISTENING**